# Technovit® EPOX

# Technovit® EPOX

- DE** Gebrauchsinformation
- GB** User Instructions
- FR** Mode d'emploi
- NL** Gebruiksaanwijzing
- ES** Modo de empleo
- IT** Istruzioni per l'uso
- RU** Инструкция по применению

# KULZER

MITSUI CHEMICALS GROUP

**KULZER GmbH**  
**Leipziger Straße 2**  
**63450 Hanau, Deutschland**  
**Tel. +49 (0) 6181/9689-2574 o. 2571**  
**Fax +49 (0) 6181/9689-2964**  
**technik.wehrheim@kulzer-dental.com**  
**www.kulzer-technik.de**

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## GB User Instructions Technovit EPOX

Technovit EPOX is a Epoxy-Resin-System, consisting of a resins component „**Technovit EPOX Resin**“ and two different hardener components „**Technovit EPOX Hardener fast**“ or „**Technovit EPOX Hardener regular**“.

Area of application is embedding of materialographic samples, especially porous materials, which have to be embedded under vacuum. Technovit EPOX characterizes itself through high transparency, strong adhesion on the sample as well as low shrinkage and UV-stability.

### Sample preparation

Clean samples carefully! The sample surface has to be free of dust, grease or other contamination. Otherwise curing can be disturbed or bubbles appear on the sample surface.

### Measuring and mixing of the components

In order to achieve an optimal result, the mixing ratio has to be adhered exactly. The ratio is two parts of „Resin“ to one part of „Hardener“ by weight.

Place „Technovit EPOX Resin“ in a mixing jar. Afterwards weight in the „Technovit EPOX Hardener“. The composition has to be mixed until the material is free of streaks (stir minimum for 60 sec.). Don't stir too heavy in order to avoid inclusion of bubbles in the material.

After mixing, let set down the material for round about 60 sec. to allow leaking the bubbles. Now pour slowly the material over the specimen, prepared in the embedding mould. Take care that there remains no air under the bottom of the specimen.

### Quantities and curing times refer to separate chart!

### Storage

Store the materials in a well ventilated place. Do not store at temperatures exceeding 25 °C. Keep away from children. Do not use after the expiry date.

### Delivery units:

Art.-No. 66040437 –  
Technovit EPOX Resin,  
1 x 1000 g

Art.-No. 66040438 –  
Technovit EPOX Hardener regular,  
1 x 500 g

Art.-No. 66040439 –  
Technovit EPOX Hardener fast,  
1 x 500 g

### Danger/Safety advice

Please follow the instructions on the product packaging and safety data sheets.  
www.kulzer-technik.com

Our technical advice, whether verbal, in writing or by way of trials, is given in good faith but without warranty, and this also applies where proprietary rights of third parties are involved. It does not release you from the obligation to test the products supplied by us as and to their suitability for the intended processes and uses. The application, use and processing of the products are beyond our control, and therefore, entirely your own responsibility. Should, in spite of this, liability be established for any damage, it will be limited to the value of the goods delivered by us and used by you. We will, of course, provide consistent quality of our products within the scope of our General Conditions of Sale and Delivery.

\* For more information please visit [www.kulzer-technik.com](http://www.kulzer-technik.com)

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### ® =marca registrada, Kulzer GmbH, Hanau
### ® = marchio registrato, Kulzer GmbH, Hanau